

COMPLETE ECMO
SOLUTION FOR ADULT,
PEDIATRIC AND
NEWBORN PATIENTS

EC \heartsuit OO



14 DAYS ECMO MODULES AND TAILORED TUBING SET

ECMO Extra Corporeal Membrane Oxygenation is an extracorporeal technique providing both cardiac and respiratory oxygen support. This technique has proved to be effective in reducing mortality of Patients suffering from potentially reversible heart failure and/or severe acute respiratory failure unresponsive to mechanical ventilation with maximal medical and pharmacological treatment⁶.

Eurosets provides Physicians with the complete solution for ECMO Patients including:

- Both PC (PHOSPHORYLCHOLINE) coated oxygenator module, and tailored tubing set, validated for 14 days use in ADULT Patients;
- PC (PHOSPHORYLCHOLINE) coated oxygenator modules, validated for 14 days, and tailored tubing sets for PEDIATRIC and NEWBORN Patients.

The A.L.ONE ECMO oxygenator family is clinically proven to be an **efficient solution for prolonged cardiorespiratory support**^{1,2,6} and is characterized by an **excellent profile in terms of biocompatibility**³ since PC coating has been shown to **remarkably reduce platelet consumption**⁴.

Moreover, being **heparin free**, the A.L.ONE ECMO oxygenator is **suitable for patients with acute HIT** (Heparin-Induced Thrombocytopenia)⁵, as in the case of heparin-induced thrombocytopenia, removal of all sources of heparin is challenging³.




TAILORED ECMO
TUBING SET
WITH 14 DAYS
VALIDATION



PHOSPHORYLCHOLINE COATING

Phosphorylcholine or PC is the predominant lipid headgroup found within the outer layer of cell membranes. PC has a natural affinity for water and binds water tightly around itself. As a result, the outer layer of the cell membrane does not promote clot formation (thrombosis).







ADULT OXYGENATOR MODULE

TECHNICAL CHARACTERISTICS

PMP fibres	
Priming volume:	225 ml
Contact surface area:	1,81 m ²
Blood flow rate:	0,4 - 7,0 l/min
Heat Exchanger surface area:	0.08 m ²
H.E. Performance Factor:	$\eta = 0.64$ (at 4 l/min)
Material:	Stainless steel
Coating:	PC phosphorylcholine
Venous Inlet:	3/8" (9.53 mm)
Arterial outlet:	3/8" (9.53 mm)

14
days







PEDIATRIC OXYGENATOR MODULE

TECHNICAL CHARACTERISTICS

PMP fibres	
Priming volume:	190 ml
Contact surface area:	1,35 m ²
Blood flow rate:	0,3 - 4.0 l/min
Heat Exchanger surface area:	0.08 m ²
H.E. Performance Factor:	$\eta = 0.64$ (at 4 l/min)
Material:	Stainless steel
Coating:	PC phosphorylcholine
Venous Inlet:	3/8" (9.53 mm)
Arterial outlet:	3/8" (9.53 mm)








14
days



NEW BORN OXYGENATOR MODULE

TECHNICAL CHARACTERISTICS

PMP fibres	
Priming volume:	90 ml
Contact surface area:	0,69 m ²
Blood flow rate:	0,2 - 1,5 l/min
Heat Exchanger surface area:	0.04 m ²
H.E. Performance Factor:	$\eta = 0.77$ (at 1,5 l/min)
Material:	Stainless steel
Coating:	PC phosphorylcholine
Venous Inlet:	1/4" (6,35 mm)
Arterial outlet:	1/4" (6,35 mm)



14
days

§ ISO 7199 condition

Bovine blood
Hb = 12 ± 1g/dL

SV02 = 65 ± 5%
PVC02 = 45 ± 5 mmHg

B.E. = 0 ± 5 mEq/L
Temp = 37 ± 1°C

Qb/Qg = 1
FIO2 = 100%



ORDERING GUIDE

CODE	DESCRIPTION	N°/PACK
EU5019	ADULT A.L.ONE ECMO 14 days (Totally Coated)	2
EU5045	PEDIATRIC A.L.ONE ECMO 14 days (Totally Coated)	2
EU5056	NEW BORN A.L.ONE ECMO 14 days (Totally Coated)	2
EU2257	HOLDER FOR A.L.ONE ECMO	1

REFERENCES

1. S.Calvi et al. **55 days ECMO support for post-cardiotomy, by-ventricular failure: a case report.** Abstract presented at 14th FECECT congress, 2011.

2. Nagi el Shouki, M.Stiller. **Case report of 19 days extra corporeal life support with a new ECMO oxygenator.** Abstract presented at 14th FECECT congress, 2011.

3. M.Pieri et al. **A new phosphorylcholine-coated polymethylpentene oxygenator for extracorporeal membrane oxygenation: a preliminary experience.** Perfusion 2012; Volume: 28 issue: 2, page(s): 132-137.

4. R. Lorusso et al. **Effects of phosphorylcholine coating on extracorporeal circulation management and postoperative outcome: a double-blind randomized study.** Interactive CardioVascular and Thoracic Surgery 8 (2009) 7–11.

5. F.Pappalardo et al. **Bioline® heparin-coated ECMO with bivalirudin anticoagulation in a patient with acute heparin-induced thrombocytopenia: the immune reaction appeared to continue unabated.** Perfusion 2009; 24: 135–137

6. M.Belliato et al. **A brief clinical case report of real-time monitoring of oxygenator performances and patient-machine interdependency during a prolonged VvECMO support.** J Clin Monit Comput. 2017 Oct;31(5):1027-1033

EUROSETS EUROPE



ITALY
Eurosets S.r.l.
Strada Statale 12, 143
41036 Medolla (MO)
Italy
Ph: +39 0535 660311
Fax: +39 0535 51248
info@eurosets.com

BENELUX
Eurosets Benelux
Winstar Park,
Rue Provinciale, 62
1301 Bierges (Wavre)
Belgium
Ph: +32 10237993
info.nl@eurosets.com
info.be@eurosets.com

FRANCE
Eurosets France Sarl
Immeuble LE FONTENOY
96, Boulevard Vivier
Merle
69423 LYON Cedex 03
Ph: +33 962191542
Fax: +33 970066573
info.fr@eurosets.com

GERMANY
Eurosets GmbH
Industriestr. 29
82194 Gröbenzell
Ph: +49 81426671227
Fax: +4981426671229
info.de@eurosets.com

EUROSETS ASIA
CHINA
Eurosets (Suzhou) Medical Device Ltd.
Suite 1308-2, The Summit (South Block)
No. 119, Suzhou Avenue West,
Suzhou Industrial Park - 215000, Greater China.
苏州工业园区苏州大道西119
号苏悦广场南楼1308-2, 215000,
更大的中国
info.cn@eurosets.com



eurosets.com

Prodotto da
Manufactured by
Fabriqué par
Hersteller
Fabricado por

Eurosets s.r.l.
Strada Statale 12, n°143
41036 Medolla (MO) Italy
Tel: +39 0535 660311
Fax: +39 0535 51248
E-mail: info@eurosets.com
eurosets.com

STERILE

0123

Latex free



Distribuito da
Distributed by
Distribué par
Vertrieb
Distribuido por

MM3EOMBE02 - © Copyright 2019 by Eurosets S.r.l. - Catalogue, products, pictures may differentiate from the actual product appearance. Specifications are subject to change without prior notice.